# TABLET FOUR-USAGE TESTER

#### **DESCRIPTION**

Tablet Four-Usage Tester is a basic equipment used for quality control of solid dosage forms. This device has four tests: dissolution, disintegration, friability, and hardness. It can be operated easily, and it has a high performance, accuracy, and sensitivity. It is widely used in pharmaceutical industry, universities, research institutes, etc.

### **TECHNICAL PARAMETERS**

MODEL	MULTITEST-4
Temperature range	20 - 45°C
Temperature stability	± 0.3°C
Speed range	20 - 200 rpm
Speed accuracy	± 2 rpm
Timing range	1 - 900 min
Timing accuracy	± 0.5 min
External size	580 x 380 x 480 mm
Weight	31 kg
Power	600 W
Power supply	AC 110V, 60Hz



Ecoshel reserves all the rights to changes in product design, if there are any design change, we will not inform in advance.

# DISINTEGRATION

Nacelles	2
Nacelle glass tubes	6
Volume of beaker	1000 mL
Nacelle up and down frequency	30-32 times/minute
Nacelle up and down range	59±2 mm



Basket units	3
Paddle units	3
Basket radial run-out	± 1.0 mm
Paddle radial run-out	± 0.5 mm
Centering	± 2.0 mm

Ecoshel reserves all the rights to changes in product design, if there are any design change, we will not inform in advance.



# FRIABILITY



Drum	1
Drum diameter	Φ 286 mm
Drum depth	39 mm
Tablet fall height	156 mm
Rotation speed	25 rpm
Rotation speed accuracy	± 1 rpm

-----

# HARDNESS

Hardness range	2 - 200 N
Probe procession range	Approx. 1 - 20 mm

#### **KEY FEATURES**

- 1. Heating, dissolution, timer, and friability functions with indicator lights.
- 2. Four LED displays are used to individual set velocity, time, and temperature and to view measurement of hardness.
- 3. Bath box is made of transparent acrylic.
- 4. Efficient temperature and circulation systems.

# **STANDARD ACCESSORIES**

Beaker

1. 2.

4.

8.

9.

- Nacelles units
- 3. Mercury thermometer
  - Bath box
- 5. Bulb syringe
- 6. Plastic hose
- 7. Brush
  - Vessels
  - Paddle units
- 10. Rotating basket units
- 11. Covers
- 12. Filter paper
- 13. Stainless steel sampling cannula
- 14. Syringe filter

Ecoshel reserves all the rights to changes in product design, if there are any design change, we will not inform in advance.