DISINTEGRATION TESTER

DESCRIPTION

2-Basket Disintegration Tester is widely used for analysis of solid dosage forms. It consists of control system, transmission system, constant temperature system, nacelle units and other accessories. It has different features as simplified operation, high accuracy and sensitivity, and stable performance.

TECHNICAL PARAMETERS

MODEL

DISITEST-2

AC 110V, 60Hz

| N | a | C | е | П | le | S |
|---|---|---|---|---|----|---|
| | | | | | | |

Nacelle glass tubes

Volume of beaker

Nacelle up and down frequency

Nacelle up and down range

Timing range

Temperature range

Temperature accuracy

Overall size

Weight

Power

Power supply

| · · · · · · · · · · · · · · · · · · · | |
|---------------------------------------|---|
| 2 | |
| 6 | |
| 1 L | |
| 30-32 times/minute | |
| | |
| 59 ± 2 mm | |
| 1-900 min | |
| 20°C to 45°C | |
| ± 0.5°C | |
| 400 x 320 x 440 mm | |
| 16 kg | |
| 600W | ! |
| | |



Ecoshel reserves all the rights to changes in product design, if there are any design change, we will not inform in advance.

STANDARD ACCESSORIES

- 1. Beakers
- 2. Nacelles units
- 3. Mercury thermometer
- 4. Bath box
- 5. Bulb syringe
- 6. Plastic hose

KEY FEATURES

- 1. Three LED displays are used to individual set nacelles and temperature.
- 2. Bath box is made of transparent acrylic.
- 3. Materials are high quality, equipment is easy-to-use.
- 4. Analysis of 6 samples for each nacelle.
- 5. Constant temperature system is composed of heater, and water pump.



Ecoshel reserves all the rights to changes in product design, if there are any design change, we will not inform in advance.