PESTICIDE RESIDUE METER

DESCRIPTION

Pesticide residue meter is portable, compact size and convenient to carry, adopts enzyme value method and shows the result of value. The pesticide residue is out of limits if 50% is positive, the higher of the value, greater of the amount of residue. It is widely used in agriculture industry, quality supervision, and agricultural training.

TECHNICAL PARAMETERS

MODEL	PESTICIDE-I
Display	LCD screen 128 x 64
Transmittance accuracy	±3%
Transmittance repeatability	±2%
Linear error	±5%
Data storage	500 records
PC interface	USB
Weight	Approx. 5.3 kg
Power supply	Four 1.5V alkaline batteries



Ecoshel reserves all the rights to changes in product design, if there are any design change, we will not inform in advance.

STANDARD ACCESSORIES



- 1. Carrying case and keys
- 2. USB cable
- 3. AC power adapter
- 4. CD-ROM
- 5. Pipette (500 to 2500 µL)
- 6. Pocket scale
- 7. Plastic test tube rack
- 8. Beaker
- 9. Glass test tubes
- 10. Glass funnel
- 11. Filter paper
- 12. Scissors
- 13. Reagents: Extract, enzyme, substrate y chromogenic agent in powder
- 14. Reagent bottle
- 15. Wash bottle
- 16. Pipette tips
- 17. Box
- 18. Extraction bottle

Ecoshel reserves all the rights to changes in product design, if there are any design change, we will not inform in advance.