PATHOLOGY WORKSTATION

DESCRIPTION

Pathological workstation is widely used in hospital pathology department, pathology laboratory, etc. Its ventilation system protects operator away from the harmful fumes produced by the formalin during sampling. The hot and cold-water system ensures it can adapt work in different climates. It is reliable quality, elegant appearance and fully functional, and equipped with high quality parts.

ECCSHEL

TECHNICAL PARAMETERS

MODEL QCT-1800

Work Surface Height

Exhaust Airflow Volume

UV lamp

LED lamp

Additional Illumination

External Size

Working Area Size

Gross Weight

Material

Consumption

Power Supply

900 mm

500-600 m³/h

20 W, emission of 253.7 nm

21 W

5 W LED gooseneck spotlight

1800 x 750 x 1950 mm

1695 x 430 x 900 mm

340 kg

Work Table: 304 stainless steel (S.S.)

Main Body: 304 S.S.

300 W

AC 110V±10%, 60Hz



Ecoshel reserves all the rights to changes in product design, if there are any design change, we will not inform in advance.

STANDARD ACCESSORIES

- 1. Adjustable spotlight
- 2. UV lamp
- 3. LED lamp
- 4. Waterproof sockets
- 5. S.S. water sink and tap
- 6. S.S. hangers
- 7. Active carbon filter
- 8. Pipe strap
- 9. Exhaust duct

KEY FEATURES

- 1. Suction inlet adopts louver-type. The upper mouthing structure ensures above blower suction, avoiding air back-streaming.
- 2. Automatic flushing device, slightly right-side inclined worktable, easy to discharge of waste and cleaning.
- 3. Base cabinet is inward-concave type, convenient for sitting operation.
- 4. Exhaust air system, high-power suction device, low noise.
- 5. Optional pulverizer can completely smash the waste during test, to avoid blocking drainage pipes.

OPTIONAL ACCESSORIES

- Pulverizer
- 2. Electrical water heater
- 3. Automatic induction hand dryer
- 4. Eye washer
- 5. Camera
- 6. Flex magnifier light
- 7. Laptop bracket

Ecoshel reserves all the rights to changes in product design, if there are any design change, we will not inform in advance.