DISSOLUTION TESTER

DESCRIPTION

6-Vessel Dissolution Tester is widely used for analysis of solid dosage forms. It consists of control system, temperature system, paddle units, rotating basket units and other accessories. It has different features as simplified operation, high accuracy and sensitivity, low noise, and stable performance.

TECHNICAL PARAMETERS

MODEL

DISOTEST-6

Basket units
Paddle units
Temperature range
Temperature stability
Speed range
Speed accuracy
Timing range
Timing accuracy
Paddle radial run-out
Basket radial run-out
Centering
Overall dimension
Weight
Power
Power supply

6
6
20°C to 45°C
± 0.5°C
20 - 200 rpm
± 2 rpm
1 - 999 min
± 0.1%
± 0.5 mm
± 1.0 mm
± 1.0 mm
960 x 320 x 480 mm
45 kg
1200W
AC 110V, 60Hz





Ecoshel reserves all the rights to changes in product design, if there are any design change, we will not inform in advance.

STANDARD ACCESSORIES

ICCSHILL)

- 1. Vessels and covers
- 2. Paddle units
- 3. Rotating basket units
- 4. Mercury thermometer
- 5. Syringe filter
- 6. Filter paper
- 7. Stainless steel sampling cannula
- 8. Bath box
- 9. Plastic hose

KEY FEATURES

- 1. Three LED displays are used to individual set velocity, time, and temperature.
- 2. Bath box is made of transparent acrylic.
- 3. Simultaneous analysis of 6 samples.
- 4. Control panel is simple and easy-to-use.
- 5. Constant temperature system is composed of heater, and water pump.

Ecoshel reserves all the rights to changes in product design, if there are any design change, we will not inform in advance.